

  
**PROF**ZrCrown



The Best  
Zirconia



PEDIATRIC  
**ZIRCONIA CROWNS**

new generation..



Made in TÜRKİYE

# About Our Zirconia



As a result of ProfTeknoloji's hundreds of model scans, more than 300 cases and 2 years of R&D work, ProfZrCrown has taken its final form.

ProfZrCrown offers a new coding technique that does not disturb the tongue. With its unique numbering system and inscription in the crown, it also eliminates the possibility of children to acquire tongue habits in the embossed part of the inscription after the application and easier tooth selection.

## WHY PROFZRCROWN?



**Innovative Wide and Inclusive Interior Design**



**Size Variety**



**Increased Durability**



**Easy Application**



**Aesthetic Design**



**Individual Coding**

# The History Of Our Company

# ProfTeknoloji

As ProfTechnology, we focus on R&D studies.

With the intraoral scanners in our clinics, the connection between the laboratory and the clinic is provided. Thanks to this connection, the patient's measurement is taken digitally and transmitted to the laboratory in a digital environment.

**2019** Company Establishment

**2020** High-end Machinery Supply

**2020  
2021** R&D Studies

**2015** Certification and Production Permit

**2015** Sales & Marketing




New  
**Product**



Our lower anterior product, which was developed with special techniques, it is now available as a set. In order to reflect your differences, you can obtain our lower anterior product as a set or refill.


*Lower*  
**Anterior**



Sizing  
**Information**

# Lower Anterior

Lower Anterior



	U1	
	U2	
	U3	
	U4	
	U5	
	U6	

Lower Anterior



	X <sub>mm</sub>	Y <sub>mm</sub>	
<b>U1</b>	<b>4.70</b>	<b>3.72</b>	<b> </b>
<b>U2</b>	<b>4.92</b>	<b>3.99</b>	<b> </b>
<b>U3</b>	<b>5.15</b>	<b>4.26</b>	<b> </b>
<b>U4</b>	<b>5.35</b>	<b>4.39</b>	<b> </b>
<b>U5</b>	<b>5.60</b>	<b>4.63</b>	<b> </b>
<b>U6</b>	<b>5.92</b>	<b>4.82</b>	<b> </b>

# Sizing Information


## Molar

**1st Upper Molar**



D2	D2
D3	D3
D4	D4
D5	D5
D6	D6
D7	D7


**1st Lower Molar**



D2	D2
D3	D3
D4	D4
D5	D5
D6	D6
D7	D7




	X <sub>mm</sub>	Y <sub>mm</sub>
D2	4.40	6.81
D3	4.60	7.05
D4	4.80	7.38
D5	4.98	7.75
D6	5.20	8.05
D7	5.50	8.45




	X <sub>mm</sub>	Y <sub>mm</sub>
D2	4.25	7.65
D3	4.55	8.02
D4	4.76	8.50
D5	5.05	8.85
D6	5.30	9.12
D7	5.74	9.60

**2nd Upper Molar**



E2	E2
E3	E3
E4	E4
E5	E5
E6	E6
E7	E7

**2nd Lower Molar**



E2	E2
E3	E3
E4	E4
E5	E5
E6	E6
E7	E7



	X <sub>mm</sub>	Y <sub>mm</sub>
E2	5.14	8.80
E3	5.30	9.11
E4	5.61	9.49
E5	5.92	9.90
E6	6.10	10.40
E7	6.47	10.96



	X <sub>mm</sub>	Y <sub>mm</sub>
E2	4.62	9.05
E3	4.80	9.36
E4	5.03	9.71
E5	5.32	10.20
E6	5.50	10.75
E7	5.75	11.20

Sizing  
Information

# Central-Lateral



**Central**

1XS	1XS	XSmall	
1S	1S		Small
1M	1M		Medium
1L	1L		Large
1XL	1XL		XLarge



**Lateral**

2XS	2XS	XSmall	
2S	2S		Small
2M	2M		Medium
2L	2L		Large
2XL	2XL		XLarge

		X <sub>mm</sub>	Y <sub>mm</sub>			X <sub>mm</sub>	Y <sub>mm</sub>	
Central		xS	5.55	6.05	Lateral	xS	5.05	5.10
	S	5.70	6.22	S		5.20	5.30	
	M	5.90	6.50	M		5.55	5.60	
	L	6.30	6.85	L		5.90	6.00	
	xL	6.55	7.22	xL		6.20	6.30	

Sizing

# Information


## Upper Canine - Lower Canine

**Upper Canine**



3XS	3XS	XSmall
3S	3S	
3M	3M	
3L	3L	
3XL	3XL	

**Lower Canine**



3XS	3XS	XSmall
3S	3S	
3M	3M	
3L	3L	
3XL	3XL	

Upper Canine



	X <sub>mm</sub>	Y <sub>mm</sub>
xS	5.75	6.32
S	6.00	6.55
M	6.25	6.90
L	6.55	7.20
xL	6.73	7.50



Lower Canine



	X <sub>mm</sub>	Y <sub>mm</sub>
xS	5.25	5.30
S	5.45	5.50
M	5.77	5.85
L	6.04	6.05
xL	6.25	6.32



New  
Case





# Box Contents

12 Pieces Central-Lateral Box

## Santral



1 Pieces



1 Pieces



1 Pieces



1 Pieces



1 Pieces



1 Pieces

## Lateral



1 Pieces



1 Pieces



1 Pieces



1 Pieces



1 Pieces



1 Pieces



Central - Lateral

36 Pieces Central-Lateral Box

## Santral



2 Pieces



4 Pieces



3 Pieces



3 Pieces



4 Pieces



2 Pieces

## Lateral



2 Pieces



4 Pieces



3 Pieces



3 Pieces



4 Pieces



2 Pieces



Central - Lateral

12 Pieces Canine Box

## Max. Canine



1 Pieces



1 Pieces



1 Pieces



1 Pieces



1 Pieces



1 Pieces

## Man. Canine



1 Pieces



1 Pieces



1 Pieces



1 Pieces



1 Pieces



1 Pieces



Max. Canine - Canine

36 Pieces Canine Box

## Max. Canine



2 Pieces



4 Pieces



3 Pieces



3 Pieces



4 Pieces



2 Pieces

## Man. Canine



2 Pieces



4 Pieces



3 Pieces



3 Pieces



4 Pieces



2 Pieces



Max. Canine - Canine

12 Pieces Lower Anterior Box

Lower Anterior

- U1**  
1 Pieces
- U2**  
1 Pieces
- U3**  
1 Pieces
- U4**  
1 Pieces
- U5**  
1 Pieces
- U6**  
1 Pieces

- U1**  
1 Pieces
- U2**  
1 Pieces
- U3**  
1 Pieces
- U4**  
1 Pieces
- U5**  
1 Pieces
- U6**  
1 Pieces



36 Pieces Lower Anterior Box

Lower Anterior

- U1**  
2 Pieces
- U2**  
3 Pieces
- U3**  
4 Pieces
- U4**  
4 Pieces
- U5**  
3 Pieces
- U6**  
2 Pieces


- U1**  
2 Pieces
- U2**  
3 Pieces
- U3**  
4 Pieces
- U4**  
4 Pieces
- U5**  
3 Pieces
- U6**  
2 Pieces






# Box Contents

84 Pieces Central-Lateral-Canine Box



136 Pieces Central-Lateral-Canine-Molar Box



48 Pieces Molar-Molar Box

Number 54						Number 64					
1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece
D7	D6	D5	D4	D3	D2	D2	D3	D4	D5	D6	D7
1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece
D7	D6	D5	D4	D3	D2	D2	D3	D4	D5	D6	D7
Number 84						Number 74					

Number 55						Number 65					
1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece
E7	E6	E5	E4	E3	E2	E2	E3	E4	E5	E6	E7
1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece	1 Piece
E7	E6	E5	E4	E3	E2	E2	E3	E4	E5	E6	E7
Number 85						Number 75					


96 Pieces Molar-Molar Box

Number 54						Number 64					
2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece
D7	D6	D5	D4	D3	D2	D2	D3	D4	D5	D6	D7
2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece
D7	D6	D5	D4	D3	D2	D2	D3	D4	D5	D6	D7
Number 84						Number 74					

Number 55						Number 65					
2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece
E7	E6	E5	E4	E3	E2	E2	E3	E4	E5	E6	E7
2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece	2 Piece
E7	E6	E5	E4	E3	E2	E2	E3	E4	E5	E6	E7
Number 85						Number 75					

# Global Market Share

Our progress in global market share continues with firm steps. Our progress in global market share continues with firm steps. We currently have a share in Turkey, which is our production location. We have market share in Dubai, Russia, Azerbaijan and Kazakhstan, respectively.



We aim to export to all countries with our ISO quality certificates. ProfZrCrown is preferred due to the features that distinguish it from other brands. We aim to increase the use of ProfZrCrown at a high rate in all countries.

We are happy with the good comments of the families of the pediatric patients as well as the physicians who use our product. Providing this service not only locally, but also globally takes us to the next level.

# Our Parent Testimonials

## Toprak's Mom

A month ago, we went to see Dr. Ozan for my 5 year old daughter's dental treatment. All our questions and our worries were over after we talked to him, and we had her treated under general anesthesia with a peace of mind. Endless thanks to Ozan who loves children very much and his job. All efforts made with love and interest. It is always very special and unforgettable for us. We also thank ProfDent employees for their good services.

15 days ago, we had my 5 year old daughter's dental treatment under general anesthesia and we were very satisfied. Ozan Baris Aydin is really good doctor. We thank him for his never-ending interest and smiling face from the beginning. Frankly, I have never seen a doctor as caring and smiling as Mr. Ozan in any private hospital before. It is a good thing that we came across him. Thanks to this process that my daughter was'nt traumatized, we cannot thank you enough.

## Ayse's Mom

## Buse's Mom

Specialist Dr. Ozan Barış AYDIN, with his smiling face and the perfection of his extensive work experience, I saw that he is a competent dentist and thank you very much, he finished 16 teeth of our daughter (6 years old) in a single session within a few hours and discharged him 1-2 hours later, his professionalism and You show that you are enough and you are indispensable with your approach that always gives information, my teacher :) without hesitation, I recommend Ozan to everyone with honor and pride.

Thanks to Profdent that our 4 year old son's dental treatment was succesfull. He went under the general anesthesio treatment, he was knocked unconscious for 3 hours and treated by our doctor Dr. Ozan. After 3 hours of hospital control, we were able to checkout after hospital. Afterwards, there was no pain. We would also like to thank our doctor Ozan Barış AYDIN.

## Aras's Mom





## Innovative Wide and Inclusive Interior Design

“ Pediatric zirconium crowns are zirconium material used in a dental treatment method for children in the field of dentistry. These types of crowns are used for the repair or replacement of children's teeth and are generally suitable for decayed, broken or extracted teeth.

[www.profzrcrown.com](http://www.profzrcrown.com)

0850 241 90 36

